

### Sens 365 Web





This user manual is intended to assist with the proper use of the Sens365 web. The Sens365 web connects to the sensor users of the CareSens Air Continuous Monitoring System (CareSens Air CGM System).

Read this user manual thoroughly and adhere to all precautions before use.

#### O Note

Thoroughly review this user manual before using the Sens365 web. Sens365 is a secondary viewer for the CareSens Air CGM System and is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice. Sens365 does not guarantee the real-time reception of glucose data of CareSens Air CGM System user.

#### i-sens

#### **Table of contents**

Be	fore us	e	6
		Document conventions	
1	Under	rstanding the Sens365 web	10
		Working Mechanism	10
2		mmended specifications for the 865 web	11
	2.1	Recommended Operating System	11
3	Using	the Sens365 web (1)	12
	3.1	Logging in	13
		Running the Sens365 web	13
		Creating an account	14
		Resetting the password	16
		Logging in	18
		Selecting the user mode	19
	3.2	Viewing the sensor user list (Care Provider Mode)	20
		Viewing the sensor user report from the list	20
		Adding favorites	
		Adding a dashboard	
	3.3	Managing groups (Care Provider Mode)	23
		Adding groups	23
		Editing group names	24
		Adding sensor users	25
		Inviting sensor users	26
		Inviting care providers	28
		Managing care providers	
		Deleting groups	31

I CODE	'n
	٠

	3.4	Report	32
		Viewing reports	32
	3.5	Settings	35
		Changing profile information	35
		Managing i-SENS accounts	36
		Resetting the password	37
		Managing connected services	38
		Deleting accounts	39
		Changing the user mode	40
		Sharing my data	41
		Accepting invites	
		Changing languages	
		Setting the threshold per diabetes type	
		Setting up AGP report criteria	47
4	•	ı the Sens365 web (2) - Sens Air Receiver	49
	Ource	Delis All Treceivel	
	4.1	Connecting the CareSens Air Receiver to the PC	
			48
	4.1	Connecting the CareSens Air Receiver to the PC	48 50
	4.1	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information	<b>48</b> <b>50</b>
	4.1 4.2	Connecting the CareSens Air Receiver to the PC Understanding the Receiver item screen	5050
	4.1 4.2	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions	<b>50</b> 5052
	4.1 4.2	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data	<b>48</b> 50505252
	4.1 4.2	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data  Send CS data	50525252
	4.1 4.2	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data  Send CS data  Changing receiver settings	<b>50</b> 5052545254
Арр	4.1 4.2 4.3	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data  Send CS data  Changing receiver settings  Customizing and resetting the receiver	<b>50</b> 505252545657
	4.1 4.2 4.3 Dendix	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data  Send CS data  Changing receiver settings  Customizing and resetting the receiver  Resetting Personal Information	48 50 52 54 56 57 59
	4.1 4.2 4.3 Dendix	Connecting the CareSens Air Receiver to the PC  Understanding the Receiver item screen  Basic Information  Exploring Receiver item functions  Uploading receiver data  Send CS data  Changing receiver settings  Customizing and resetting the receiver  Resetting Personal Information  A Frequently Asked Questions	48 50 52 54 56 57 59

#### i-sens

#### **Before use**

- This product only makes use of personal information which the user has agreed to allow to be collected.
- Contact the manufacturer or visit the website for details of this product.
- Check whether the device used to run the software is turned on and connected to a wired/wireless network.

#### O Note

- All standards and regulations mentioned are in effect as of the date this
  manual was issued.
- i-SENS, Inc. has carefully prepared the information in this manual to be as
  accurate as possible. However, i-SENS is not responsible for any errors
  or omissions contained in the manual. i-SENS may change the product
  described in this manual or any related software applications without
  notice in order to enhance the product reliability, features, or design.
- This document is protected by copyright. It is strictly prohibited to copy or alter this manual without prior consent from i-SENS.



#### **Document conventions**

#### **Definitions of symbols**

These symbols not only provide additional information on the product and product use, but also on safety.

Symbol	Description
CE	CE Mark
EU REP	Authorised representative in the European Community/ European Union
MD	Medical device
[]i	Consult instructions for use or consult electronic instructions for use
***	Manufacturer
#	Model number
UDI	Unique device identifier



#### **Safety Information**

#### Indication for use

Sens365 is indicated for continuous monitoring of blood glucose levels via measurement of glucose in the interstitial fluid in persons with diabetes mellitus aged 18 years and older. Sens365 is intended for use by patients, care providers, and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis, and evaluation of historical glucose data from compatible glucose monitoring devices to support effective diabetes management. Sens365 is available as cloud-based web software (Sens365 Web) and as a smart device application (Sens365 App).

Sens365 Web is not intended to provide clinical decisions or to be used as a substitute for professional healthcare advice. Patients should consult a healthcare professional before making any medical interpretation and therapy adjustments based on the information provided by the software. Healthcare professionals should use the information in conjunction with other relevant clinical information available to them.

Sens365 App is a secondary viewer that allows home users and healthcare professionals to remotely view glucose data from users of compatible glucose monitoring devices. The software is intended solely to provide a secondary level of awareness and should not be expected to communicate and display real-time glucose information.

#### **Important User Information**

- You must create an i-SENS account to use Sens365.
- Sens365 is a secondary viewer for the CareSens Air CGM System and is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice. Sens365 does not guarantee the realtime reception of glucose data of CareSens Air CGM System user.
- The warnings and cautions must be read and understood before using Sens365, and they must be strictly adhered to. Sens365 must be used according to the instructions in the user manual.



#### **Warnings**

- Users of CareSens Air CGM System and their family members and guardians must not make treatment decisions based on the Sens365 glucose data. In addition, Sens365 does not offer any medical advice and must be not used for that purpose. You must consult a medical professional to make a treatment decision.
- Medical institutions and healthcare professionals should use data provided by Sens365 in conjunction with other clinical information available to them for treatment decisions or analysis.

#### **Precautions**

- CareSens Air CGM System user must turn on the share feature in order to send glucose data to his or her care providers (family members, guardians, medical institutions, medical personnel). If the care providers want to view the shared content, they must download Sens365 and connect with CareSens Air CGM System user.
- The care providers' smart devices must be connected to the Internet in order to view the information shared by the CareSens Air CGM System user through Sens365.
- If the smart device of the CareSens Air CGM System user or care provider (family members, guardians, medical institutions, medical personnel) is turned off, in airplane mode or not connected to the Internet, the glucose data may not be sent or received.
- Care providers (family members, guardians, medical institutions, medical personnel) can co-manage the CareSens Air CGM System with the user, however, the sensor user must always be the subject of active glucose management.
- If care providers (family member, guardian, medical institution, medical personnel) terminate the connection, the content sent by CareSens Air CGM System user can no longer be viewed.

#### i-sens

#### 1 Understanding the Sens365 web

#### **Working Mechanism**

This product is a secondary display software that collects and saves processed glucose data with i-SENS' CareSens Air CGM System. The system is capable of sending data and displaying information on Android and iOS smart devices. The glucose data collected and processed by the i-SENS' CareSens Air CGM System can be securely transmitted to care providers (including family members, guardians, medical institutions, and medical personnel) using the Secure Sockets Layer (SSL) encryption protocol.

# 2 Recommended specifications for the Sens365 web

The following operating systems have been verified as compatible by their manufacturers.

## 2.1 Recommended Operating System

#### O Note

Make sure to use an operating system that meets the recommended specifications. Proper functioning of the Sens365 web is not guaranteed for operating systems which have not been verified as being compatible.

Recommended Browser	Resolution
Google Chrome, Mozilla Firefox, Microsoft Edge, Apple Safari	1280 x 800 px or higher

OS (Operating System)	Version
Windows	10 or later
Mac OS	11 or later

#### i-sens

## 3 Using the Sens365 web (1)

You can run the Sens365 web and monitor the glucose readings measured by the connected continuous glucose monitoring system. In order to log in to the Sens365 web, you must create an account.

Using various features offered by the Sens365 web enhances the efficiency of diabetes management.

This section will help you to:


- 1. Register for and log in to the Sens365 web.
- 2. View the sensor user list.
- 3. Manage groups, sensor users, and care providers.
- 4. View the Reports screen and glucose reports.
- 5. Edit basic information and manage i-SENS accounts by changing settings.



#### 3.1 Logging in

#### Running the Sens365 web

- Open your browser and enter the Sens365 web url. URL: https://sens365.com
- 2 The Sens365 web start screen is displayed.



#### **Creating an account**

If you do not have a Sens365 account, go to 'Create an account' screen of the login screen to sign up.

- 1 Click **Get Started** on the Sens365 web start screen.
- 2 Click Create an account when the account login screen is displayed.
  - Even if you choose Social Login, you can sign up with an i-SENS account.
  - To sign up with Social Login, enter your country of residence, and user information, and agree to the required terms and conditions.



#### O Note

Change the language on the top right as needed.

3 Verify your email address.





4 After accepting the required terms and conditions, click **Ok**.



5 Enter your password and fill out all your information, then click Ok. Your account registration will be completed successfully.



#### Resetting the password

If you forget your password or need to change your password for security reasons, you can reset it.

- 1 Click **Get Started** on the Sens365 web start screen.
- 2 Click Forgot your password? when account login screen is displayed.



3 Verify your email address.





4 Set your new password.



5 Log in to your account with the new password.

#### Logging in

If you have a Sens365 account, log in to Sens365 on the login screen.

- Click **Get Started** on the Sens365 webstie start screen.
- 2 Enter your registered email and password when the account login screen is displayed.
  - If you have signed up with Social Login, Keep using the same account that was selected in Social Login.



#### Logging out

Click Log Out after selecting the profile on the top right of the screen to log out.



You will be redirected to the Sens365 start screen if you logged out properly.



#### Selecting the user mode

The user mode selection screen will be displayed when you log in to Sens365 for the first time.

- CareSens Air users should select the Sensor User mode.
- Care providers of CareSens Air should select the Care Provider mode.
- The user mode selection screen will be displayed when you log in to Sens365 web for the first time.



- Select Sensor User or Care Provider according to the use purpose.
- 3 For a sensor user, 'Report' appears on the start screen.



For a care provider, the sensor user list is displayed on the start screen.



## 3.2 Viewing the sensor user list (Care Provider Mode)

#### Viewing the sensor user report from the list

1 If you select a group on the left, the sensor user list of the selected group will be displayed.



2 If you select a sensor user in the sensor user list, you can view the report on the selected sensor user.





#### **Adding favorites**

You can add frequently searched sensor users of the sensor user list to the favorites list.

1 Click ☆ to the left of the desired sensor user from the sensor user list. The ☆ changes to yellow when added to favorites.



Check the added sensor user in the favorites list.



3 Once a sensor user has been added to the favorites list, clicking on ★ again removes them from favorites.

#### Adding a dashboard

Sens365 offers a dashboard feature to be able to quickly view the reports of multiple sensor users.

1 The dashboard can be viewed by clicking on Dashboard on the top right of the screen.



2 Click Add to Dashboard of the desired user from the sensor user list to add the sensor user to the dashboard.



3 Open 'Dashboard' and check whether the selected sensor user has been added.



- 4 Click of the sensor user to delete the sensor user from 'Dashboard'.
- 5 To clear 'Dashboard', click i of the group and Empty Dashboard.

#### Note

Click  $^{\circ}$  to pin sensor users in 'Dashborad'. Sensor users pinned in 'Dashborad' are not deleted from the list.



## 3.3 Managing groups (Care Provider Mode)

#### **Adding groups**

1 Click: in the 'My Group' on the left and click Add Group.



2 Enter the name of the group to add.



3 A group is created with the entered name.




#### **Editing group names**

1 To edit the group name, click on : of the group you would like to edit and click **Group Settings**.



2 Enter a new group name and click Save.



3 Check the edited group name.





#### **Adding sensor users**

1 To add a new sensor user to the group list, click Add on the top right of the list.



2 Enter the share code of the sensor user to add and click **Ok**.




#### O Note

The share code can be created in the CareSens Air app used by the sensor user. For more details, check the CareSens Air user manual.

3 If the entered share code is valid, the sensor user is added to the list.

#### **Inviting sensor users**

To invite a sensor user to the group, click : of the group and click Invite Sensor User.



2 Enter the group name to display to the invited sensor user and click Next.




#### O Note

The entered group name will be visible to the invited sensor user.



3 Check the created share code and click **Email Invitation**.



4 Enter the email address of the sensor user to invite and click Send. You can send a maximum of 10 emails at once.




5 An invite email is sent to the provided email address.




#### **Inviting care providers**

1 To invite care providers to the group, click : of the group and click Invite Care Provider.




2 Enter the email address of the care provider to invite and click **Send**.





#### Managing care providers

To manage the care providers of the group, click of the group and click Manage Care Providers.



2 Select the group to manage from the groups displayed in the group care provider management list.



The care providers that belong to the selected group are listed.





4 To change the permissions of the care providers of the group, click ▼ of the care provider permissions and select the permission to change.



#### O Note

There are 'Admin' and 'Coordinator' permissions.

To delete a care provider from the group, click ... on the right of the care provider to delete and click Remove.






#### **Deleting groups**

#### O Note

Groups can be deleted if there are no registered sensor users in the group.

1 To delete a group, click: of the group to delete on the left menu, and click **Delete Group**.



2 Click **Ok** when the delete confirmation pop-up appears. The selected group will be deleted.





#### 3.4 Report

#### **Viewing reports**

- 1 Follow these steps to check the report.
  - Sensor User mode: Click 🔀 Report on the top left menu.
  - Care Provider mode: Click the row of the desired user from the sensor user list.
- 2 The 'Report' screen is displayed.



#### O Note

- In Sensor User mode, the 'Overview' screen is displayed first.
- In Care Provider mode, the 'AGP' screen is displayed first.



3 Each report is composed of the following content.

Report	Content
Overview	The following information for the selected period is provided.
	Time in Ranges
	Average Glucose
	Glucose Management Indicator (GMI)
	Glucose Variability
	Standard Deviation
	Daily Average Glucose
	Insulin
	Daily Time in Ranges
	Time CGM Active
	Device Information
	• Sensor
	Min/Max Glucose by Events (fingerstick)
	Ketones
	• GKI
	Average Meal per Day
	Exercise
AGP	The standard AGP for the selected period is provided.
Comparison	A table and chart is provided to compare glucose levels per period.
	Glucose within Target Comparison
	Comparison (Glucose Variability, Standard
	Deviation, Glucose Management Indicator (GMI), Average Glucose)
	Patterns Comparison
Patterns	The daily charts and event information on areas where low glucose, time in ranges, and high glucose have occurred during the selected period are displayed.
Trends	The daily charts and weekly charts for the selected period are provided.



Report	Content
Statistics	The following information for the selected period is provided.
	Time in Ranges
	Average Glucose
	Average Glucose Variability
	Average Time CGM Active
	Hourly Time in Ranges
	Hourly Average Glucose Value
	Hourly Average Glucose Variability
	Hourly Average Time CGM Active
Events	User events entered during the selected period are provided in a table.

- 4 To save the report, click 

  on the top right.
  You can select the report and save it as a PDF file.
- 5 To print the report, click on the top right. You can select the report and print it.
- To send the report by email, click 

  on the top right. You can send the report PDF file to a specific email.



#### 3.5 Settings

#### **Changing profile information**

To change basic information, click the profile at the top right or User Settings on the left.



- The basic information that can be edited is as follows:
  - Name
  - · Date of Birth
  - Gender
  - Phone
  - Diabetes Type
- 3 When information is changed, the changed information is saved automatically.


#### Managing i-SENS accounts

1 To manage your account, click **Go to i-SENS Account Settings** on the personal information screen.



2 Click Edit Profile and check the information entered upon creating the account.







#### Resetting the password

To change your account password, click Change Password on the account settings screen.



2 Enter the current password then enter the new password.



#### Managing connected services

1 To manage connected services, click **Manage Connection Services** on the account settings screen.



You can check the list of services in use with your current account.





#### **Deleting accounts**

1 To delete your account, click Delete Account on the account settings screen.



When deleting your account, check the notice and click I agree and Ok to complete the process.




#### Changing the user mode

1 To change the use mode, click **Switch User Mode** on the personal information screen.



When the Switch User Mode pop-up is displayed, select the user mode to change to and click **Ok**.



3 The screen layout will be changed according to the selected user mode.



#### **Sharing my data**

#### O Note

You can enter the share code for glucose data sharing in 'Sensor User' mode.

To share your glucose data with a new care provider, click **Connections**. A list of care providers who receive your shared data is displayed.



- To enter the share code, click Connect with Care Provider on the top right.
- 3 Click Ok after entering the share code of the care provider to receive the shared data.



4 If the entered share code is valid, the group is added to the care providers list.




To stop sharing with an existing care provider, click **Disconnect** on the right of the care provider to disconnect.





6 Click **Ok** when the stop sharing confirmation pop-up appears. Sharing will be stopped with the selected care provider.





#### **Accepting invites**

1 If you are invited to a group as a health manager, the invitation information is displayed in 'Invitation List'.



- 2 To accept an invite, click Accept.
- 3 Click **Ok** when the accept confirmation pop-up appears.



4 The invited group is displayed on the left.





- 5 To reject an invite, click Reject.
- 6 Click **Ok** when the refuse confirmation pop-up appears. The invited group is not added to the menu.





#### **Changing languages**

- 1 To change the language, click at the bottom right of the screen.
- 2 Select the desired language on the language selection screen.



#### Setting the threshold per diabetes type

- To set the threshold for setting the target diabetes glucose range, click Target Range Settings at the top of 'Report' screen.
- 2 Set the target glucose per area on the 'Target Range Settings' screen.



#### Note

- Clicking on 'Reset' resets to the diabetes type threshold.
- Changing the target glucose range settings applies to the Time in Ranges of the sensor user. (AGP, Comparison excluded)
- 3 If a day/night application time change is needed, the Day Start Time and Night Start Time can be set up in the Day & Night Time section.



#### Setting up AGP report criteria

- 1 To set the custom range of the AGP report, click **AGP** on the 'Report' screen.
- 2 If you set the report criteria at the top of the screen to user settings, you can change the low or high criteria.



3 When you click Apply, the AGP report on the screen will be displayed based on the set value.

#### i-sens

# 4 Using the Sens365 web (2) - CareSens Air Receiver

# 4.1 Connecting the CareSens Air Receiver to the PC

CareSens Air Receiver users: follow the steps below to connect the receiver to the PC.

- 1 To connect Sens365 to the CareSens Air Receiver, go to the CareSens Air Receiver menu.
  - If you are not logged in, click Uploader Software at the bottom of the Sens365 web start page.



If you are logged in, after signing in, click Uploader Software at the bottom
of the website.






2 If Uploader Software is not installed on your PC, follow the on-screen instructions to install the version for your OS. When installation is complete, click Refresh on the right.



3 Connect the CareSens Air Receiver to the PC. Follow the on-screen instructions to complete the connection. When the connection is complete, click Ok at the lower right.





# 4.2 Understanding the Receiver item screen

#### **Basic Information**

The table below describes the basic information shown when the CareSens Air Receiver is connected to the PC and is communicating properly with the Sens365 web service.



No.	Name	Description		
0	Connected receiver information	Displays basic information such as product name, serial number, model name, and software version.		
0	Send CS data	Sends log data from the connected receiver.		
0	Sync Time	Synchronizes the CareSens Air Receiver time with the PC time.		
4	Reset Receiver	Resets the connected CareSens Air Receiver.		
6	Software Update	<ul> <li>Updates the CareSens Air Receiver software to the latest version.</li> <li>Active only when the receiver software is not up to date.</li> </ul>		
0	Personal Information	Displays the designated user information.		
0	Select User	Designates the user who uses the CareSens Air Receiver.		



No.	Name	Description	
3	Manage Data	Uploads data stored on the CareSens Air Receiver to the Sens365 web service.	
9	Settings	Allows to view and edit basic settings (language, date format, time format, etc.) and notification settings.	

#### O Note

The ' $\Theta$  Sync Time', ' $\Theta$  Reset Receiver' and ' $\Theta$  Software Update' functions are disabled and cannot be selected while the sensor is connected. These functions are enabled only when the sensor is not connected.



# 4.3 Exploring Receiver item functions

#### **Uploading receiver data**

- 1 To upload the latest data from the CareSens Air Receiver, select the Receiver on the left
- 2 Click Select User in Personal Information.



- 3 Check the **Personal information** and click the **Manage Data** tab.
  - In the Manage Data tab, the number of data records stored on the CareSens Air Receiver is displayed.



4 Click **Upload** at the bottom of the screen to upload the receiver data to the Sens365 web service.




#### O Note

The Manage Data tab is visible only to logged-in users.

#### i-sens

#### Send CS data

- To upload CS data from the CareSens Air Receiver, select the Receiver on the left.
- 2 Click Select User in Personal Information. (If a user was previously assigned, it will already be assigned.)
- 3 Check the Personal Information and click Send CS data.




In the Contact Us pop-up, write your inquiry and click Send CS data at the bottom. (To send only CS data without an inquiry, leave the field blank and click Send CS data.)






5 Click **Ok** to send the CS data.



#### **Changing receiver settings**

- 1 To change settings of the connected receiver, select the **Receiver** on the left.
- 2 Select the basic and notification setting items you want to change, then click Apply.



3 When information is changed, it is saved automatically.



#### **Customizing and resetting the receiver**

You can assign a specific user to the receiver for management.

Once a user is assigned, log data and uploads are performed under that account.

#### O Note


If the user assigned to the receiver does not match the Sens365 logged-in user or lacks permission for the assigned user, certain functions are restricted.

#### Selecting a user

Depending on the 'User Mode' selected in Sens365, the method for uploading receiver data will differ for the sensor user.

For user-mode settings, see 'Selecting the user mode' or 'Changing the user mode'.

- To assign a receiver user, click the **Receiver** on the left.
- 2 In Personal Information, click Select User.






3 For sensor users, designate the currently signed-in Sens365 account as the receiver user and click **Select**.

For healthcare professionals, designate one of the data-sharing users as the receiver user and click **Select**.

• In the 'General' group, only users with 'Selectable' status can be selected.



4 The receiver user is assigned.



#### **Resetting Personal Information**

- 1 To assign a receiver user, click the receiver on the left.
- 2 After checking Personal Information, click Reset User Information. The Personal Information will be reset.



#### i-sens

## Appendix A Frequently Asked Questions

This chapter presents situations that may occur while using Sens365 and how to deal with them. If any situation which is not presented in this chapter occurs, or if you experience an issue that you are unable to resolve on your own, contact the authorized distributor.

This section will help you to:

- Identify the causes of problems that occur while using the Sens365 web.
- Resolve problems that occur while using the Sens365 web.

#### **Dashboard**

#### Care Provider

#### How do I search for sensor users?

You can search for sensor users by names, notes, and tags included in the list.

#### How can I manage or monitor a specific sensor user?

By clicking 🚖 of a user in the sensor user list, you can set them as a favorite. Then, you can manage them separately in the favorites tab.

#### How can I filter the list of sensor users?

The sensor user list can be viewed by applying the filtering or sorting options provided in the table header of the sensor user list.

#### What is the list of sensor users?

It is a list of sensor users connected for management in the respective group. You can check the real-time glucose information of sensor users in the group.

#### **General**

#### Sensor User

#### How do I upload my glucose data?

Glucose readings in the CareSens Air app will be automatically uploaded. You can also manually upload data in the CareSens Air app > **Settings** > **Upload Now**.



However, you can only upload data when you have an internet connection.

#### What is the dormant account policy for Sens365?

If a user doesn't use Sens365 for a year, their account will be classified as dormant and their personal information will be securely stored separately. This stored information will be retained for a year and may be discarded if there's no legal requirement to keep it. i-SENS will notify dormant users 30 days before their account becomes dormant and provide information about the stored data. If you don't want your account to become dormant, simply log in before the dormancy period. Even if your account becomes dormant, as long as your personal data hasn't been discarded, you can log in, obtain user consent, restore the dormant account, and continue using the service.

### How can I ensure that the data stored in Sens365 is securely protected?

Personal and sensitive information stored in Sens365 is encrypted. During data transmission, the SSL (Secure Sockets Layer) encryption protocol ensures the confidentiality and integrity of the data, protecting it from eavesdropping and alterations.

#### **Care Provider**

#### What is the dormant account policy for Sens365?

If a user doesn't use Sens365 for a year, their account will be classified as dormant and their personal information will be securely stored separately. This stored information will be retained for a year and may be discarded if there's no legal requirement to keep it. i-SENS will notify dormant users 30 days before their account becomes dormant and provide information about the stored data. If you don't want your account to become dormant, simply log in before the dormancy period. Even if your account becomes dormant, as long as your personal data hasn't been discarded, you can log in, obtain user consent, restore the dormant account, and continue using the service.

#### Is the data of the added sensor user in the group secure?

The data is securely transmitted during the transmission through the SSL (Secure Sockets Layer) encryption protocol, and personal and sensitive information is stored in an encrypted state. Only care providers belonging to the group can access the data of the added sensor users in the group.

#### i-sens

#### Group

#### **Care Provider**

#### How can I join someone else's group?

The care provider of the group you wish to get invited to can click next to that group and click **Invite Care Provider** to send an invitation via email. By clicking the invitation email link, you can accept the invitation and join the group. Additionally, you can accept the invitation from the invitation lists in **Settings** > **Profile**. However, the inviting email must match the care provider's Sens365 login account.

#### How many sensor users or care providers can be in one group?

There is no limit to the number of sensor users and care providers that can be added to a group.

#### Who can I invite to the group?

You can invite sensor users to the group to manage their glucose levels, and you can also invite other care providers to co-manage the sensor users of the group.

#### What is the group sharing code and how can I check it?

The group sharing code is a unique code generated for each group to invite and manage sensor users. You can check the sharing code of the group by clicking in the group menu and selecting **Invite Sensor User**. Please note that each group has a unique sharing code, and once assigned, the sharing code does not change. If you want to change it, click C next to the sharing code.

#### What features are available in the group?

Care providers create groups for organizational and management purposes to effectively manage sensor users. When you first join Sens365, a default group is automatically created under your name. You can change the group name as you wish, create subgroups, and manage sensor users. Additionally, you can invite or add sensor users within the group to check their reports, and you can also invite care providers to the group to co-manage sensor users.



#### Mode

#### Sensor User

#### How do I set up Sens365 in sensor user or care provider mode?

After signing up for Sens365, a mode selection screen will appear during your first login. Select either the sensor user or care provider mode, and Sens365 will operate in that mode. You can later change the user mode in Profile section under the Settings menu.

#### Care Provider

#### How can I set up as a sensor user or care provider?

After signing up for Sens365, a mode selection screen will appear during your first log in. From this screen, you can choose the sensor user or care provider mode, and Sens365 will run in the selected mode. Later, you can change the user mode under the Settings menu in the Profile section.

#### Reports

#### **Sensor User**

#### How can I share reports?

Reports can be shared through a connection with a care provider in two ways:

- Enter the sharing code provided by your care provider at Settings > Connections > Connect New Sensor User.
- Share the sharing code generated in the CareSens Air app (CareSens Air app > Settings > Care Providers > Connect > Generate Sharing Code) with the care provider.

#### How do I check the reports?

#### **Care Provider**

#### How do I download the raw glucose data of the sensor user?

You can download the raw glucose data in CSV format by clicking **(x)** in the top right corner of the sensor user's report screen.



#### How do I compare the glucose data of the sensor user?

Glucose data can be compared by setting two periods in the Comparison tab of the sensor user's report screen.

#### How can I check the report?

When you click a sensor user within a group, you can view their glucose report. If you want to print the report, click on  $\displaylim$ , or  $\ensuremath{\boxtimes}$  in the top right corner of the report screen to view the report in PDF format.

#### **Data Sharing**

#### **Sensor User**

#### What should I do if I can't connect with a care provider?

Please verify if the sharing code received from the care provider is accurate.

#### How do I stop sharing my data?

To stop sharing your data with a care provider, go to Settings > Connections and choose to disconnect.

#### Care Provider

### How can a care provider view the glucose data uploaded by the sensor user?

The care provider must be connected with the sensor user through their affiliated group to view the glucose data uploaded by the sensor user. There are two methods for this:

- The first method is to click in the group menu where you manage the sensor users. Then sclick Invite Sensor User. Next, directly provide or email the sharing code of that group to the sensor user. The sensor user then enters the group's sharing code (CareSens Air app > Settings > Care Providers > Connect > Enter Sharing Code). Please note that each group has a unique sharing code, and once assigned, the sharing code does not change. If you want to change it, click next to the sharing code.
- The second method is to click of the group managing the sensor users, click
  Add Sensor User. and then enter the sharing code created in the CareSens Air
  app (CareSens Air app > Settings > Care Providers > Connect > Generate
  Sharing Code) to add the sensor user.

# Appendix B Technical Information

#### **B.1** Performance

Web Browser: HS-FE-01

Product name	Sens365	
Compatible medical devices	CareSens Air Continuous Glucose Monitoring System	
Communication environment	Communication protocol: HTTPS	
Operating environment	<ul> <li>Recommended browser: Google Chrome, Mozilla Firefox, Microsoft Edge, Apple Safari</li> <li>OS(Operating System) Version: Windows 10 or later, Mac OS 11 or later</li> <li>Resolution: 1280 x 800 px or higher</li> </ul>	
Transmitted and reveiced data	<ul> <li>Continuous glucose values of a compatible medical device user, user input data (meals, insulin, exercise, etc.)</li> <li>User data statistics and reports including the glucose of a compatible medical device user</li> </ul>	
Key performance	<ul> <li>Care providers (family members, guardians, medical institutions) who use the app can receive shared data by connecting with the user web of a compatible medical device user.</li> <li>Displays statistics for a compatible medical device user's data and saves reports.</li> <li>Displays the continuous glucose values of a compatible medical device user and user input data.</li> </ul>	
Software version	Version: 1.1.x	



### **B.2 Cybersecurity**

#### Cautions regarding cyber threats

- Install a virus protection or antivirus program on your smart device, whether a Windows, Mac OS or tablet device, to prevent malicious programs from accessing the device's information.
- If you are unable to receive the glucose data from the sensor users despite
  entering the share code, make sure the six-digit share code was entered
  correctly. If the same issue persists despite entering the correct share code,
  contact Customer Service via the product website, caresensair.com.
- If you are logged out due to another device logging in, change your password if it was not you who attempted to log in.
- If you enter the wrong password more than five times, reset your password.
- If you have not used your account for over 1 year, the account becomes inactive. If your account has become inactive, you must reactivate it.
- If a cybersecurity incident occurs related to Sens365, contact the nearest authorised distributor and the competent authority.

#### Safety regarding cyber threats

- A token is issued to the user once his or her credentials have been verified as valid.
- For each request, the validity of the token is verified, and then the authentication is processed.
- The SSL communication method is used to send data.
- AES encoding is used to encrypt personal information data and the data is saved to a DB.
- If you enter the wrong password more than five times, you must reset your password.
- If you have not used your account for over 1 year, the account becomes inactive. If your account has become inactive, you must reactivate it.

### i.sens



i-SENS, Inc.

43, Banpo-daero 28-gil Seocho-gu, Seoul 06646 Republic of Korea CareSensAir.com

© 2025 i-SENS, Inc. All rights reserved. PCGAA-000121 REV1 2025-11



Medical Technology Promedt Consulting GmbH

Ernst-Heckel-Straße 7 66386 St. Ingbert, Germany



